

# Calibed

Warming bed



Calibed warming bed provides the ideal and constant temperature that newborn babies need to feel secure. This protective heat also offers feather-soft comfort for baby. Calibed provides child's skin temperature monitoring. The high and low alarms prevent any risk of hypothermia or hyperthermia.

# Calibed - Ref.1091

Warming bed

## SPECIFICATIONS

**Mattress temperature regulation** from 28 to 38.5°C  
(with secure access under 35°C and above 37°C)

**Uniform warming** over the entire bedding area

**Skin temperature monitoring**

**Sound and visual alarms.** Advanced, offset management system for alarms (Quietude mode). 3 sets of visual and sound alarms have been defined depending on the level of emergency for the comfort of newborn baby and the medical team.

**Heating power :** 75 watts (absorbed power : 170 VA)

**Power supply :**

- Control box : 220-240V / 50-60Hz
- Mattress : 24VAC security low voltage

**Electrical Classification :** class I

**Class IIb medical device**

## DESCRIPTION

**Bedding built-in control box** with LCD screen for a good visibility of data and alarms

**Bedding :** 720x510mm surrounded by 4 transparent protection shields around the baby.

**Very comfortable polyurethane gel mattress** (latex- and silicone-free).

**Tilt positions :** +17/-10° continuous

**Fold down side panels** for an easy access to the baby. There are fully removable for washing by soaking.

**4 tubing ports and 4 orifices for drains** to pass through at low point.

**Tent** (optional) is intended to decrease external nuisances which may disturb the new-born (light, noises...). It can be positioned to isolate partially or totally the patient from light to recreate an alternation between day and night even to respect and promote the circadian rhythm. Removable textile for cleaning.

**Floor/bedding height :** 870/1070 mm

**4 swivelling castors** diam. 125mm (2 with brake)

**General dimensions (LxWxH) :** 600 x 1070 x 1089 mm

**Delivered with one skin probe**



## OPTIONS & ACCESSORIES

### Tent (ref.1092)

Designed to decrease the surrounding nuisances (light, noises etc.)

In cotton (OEKO TEX STD 100 CLASS I)

With flexible support to quickly adapt the position of the tent.

### Storage space (ref.4339)

2 drawers (with 1 pivoting) and a rear shelf.