



Micro-granules changing mattress

Ref. 3729

For newborn comfort

Carewave by Médipréma micro-granules changing mat fits with all changing tables. It provides very good conditions of comfort to the baby. It is an alternative to the classic foam changing mat whose contact may be cold and hard. The micro-granules mattress offers a soft contact. Nestled against the cushion, the baby feels more secure and is more quiet during cares and changes.



Easy to use and secure

The Carewave by Médipréma changing mat is light and easy to install. Made of polystyrene micro-granules in a polyurethane coating, it is easy to clean by sponging.

Specifications:

Dimensions	57 x 55 cm
Material inside the cushion	Micro-Granules in polystyrene
Matière enveloppe	Bi-elastic polyurethane slipcover
Cleaning	by sponging Water resistant Double seam
Packing	Per unit

A Pharmaouest product, distributed by: