



Foam & microbeads mattress

Special for transport

Ref : 5044

Support and comfort for newborn babies

By combining both a comfortable memory of shape foam mattress with polystyrene micro-beads acknowledged for their maintaining properties, the mattress *Carewave by Médipréma* brings an effective solution to the problems of safety and comfort during transport.



- > The memory of shape mattress is made of visco-elastic which moulds the shape of the child's body. Its soft comfort enables to reduce inconvenience caused by jolts during transport.
- > The additional micro-beads filling enables to maintain the newborn body. Better retained the baby feels more secure and more quiet during transport.

Specially sized to be used in Médipréma Nite transport incubator, the foam & micro-beads transport mattress, placed directly on the mattress, provides ideal comfort conditions to receive the baby.

The foam and micro-beads transport mattress can be attached to the bedding area thanks to the Magic Strips straps. The Magic Strips is a limp micro-fiber velcro strap, soft to the touch, very easy to put in place as well as to take off.



The foam and micro-granules transport mattress is light and easy to put in place. It can be taken off as easily in order to quickly install the baby on the hard bedding for special medical procedures.

Specifications :

Materials	<i>Mattress filling:</i> Visco-elastic polyurethane foam mattress + polystyrene micro-beads <i>Cover :</i> Bi-elastic polyurethane cover, suited to skin contact
Dimensions (LxW)	61 x 38 cm
Cleaning	Easy to clean and disinfect by sponging Water resistant thanks to double seam
Packaging	Per unit