



Microbeads positioning cushions

For newborn and infant positioning



Security and posture prevention

The Carewave by Médipréma cushions protect newborns from skin lesions thanks to their anti-bedsore function. They limit the risk of osteo-articular distortion and muscular retraction in positioning baby's members in a right axis. Composed with polystyrene microbeads, ultra flexible positioning cushions prop up newborns without press them.

Comfort and freedom of movement

Médipréma items reduce newborn's pain and stress, improving its general clinical state. By allowing voluntary movements and preventing from primary reflexes, they favour motor development. Cushions also facilitate intervention of hospital staff and contact with baby's parents.



Two models available



"Demi-lune"

The "Demi-lune" cushion allows placing easily and quickly the infant in axis lining up head, spinal vertebral and legs.

The risk of candelabra arm and the "squashed frog position" are dismissed. In lateral, ventral or dorsal decubitus, the baby is maintained with security.

"Flexitube"

The "Flexi-tube" cushion offers freedom of use in permitting to mould the cushion depending on medical cases and needs. It can also warp around baby's body and, that way, reproduce maternal cocoon.

Also available: washable jersey stitch protection cover



Microbeads positioning cushions:

- Polystyrene micro-granules
- Bi-elastic polyurethane and polyamide cover
- Water resistant
- Double seam
- Maximum weight accepted: 25 kg
- Cleaning by sponging

Benefits:

- Reduce the risk of osteo-articular distortion and muscular retraction
- High comfort and freedom of movement
- Easy and quick to use
- Resistant to time and regular uses
- Easy cleaning and disinfection

Demi-lunes references:

Patient	Demi-lune		Protection cover	
	Ordered reference	Maker reference	Ordered reference	Maker reference
Premature baby ≤ 1 kg (99 cm)	3701	NEG0699	3702	NEG0702
Premature baby 1 à 2 kg (124 cm)	3703	NEG0700	3704	NEG0703
Newborn (136 cm)	3705	NEG0701	3706	NEG0704
Infant (174 cm)	3707	NEG0705	3708	NEG0706



Flexitubes references:

Patient	Flexitube		Protection cover	
	Ordered reference	Maker reference	Ordered reference	Maker reference
Small (140 x 8 cm)	3709	NEG0707	3710	NEG0710
Medium (170 x 11 cm)	3711	NEG0708	3712	NEG0711
Large (200 x 15 cm)	3713	NEG0709	3714	NEG0712



Pharmaouest products, distributed by:



ZA Node Park Touraine - 470 rue Gilles de Gennes - BP 8 - 37310 TAUXIGNY - FRANCE
Tél +33 (0)2 47 28 47 87 - Fax +33 (0)2 47 27 35 85 - www.mediprema.com

