

Micro-granules cushions for weighing-scale

Ref. 3728

To facilitate the new-born weighing

- Less stress for the newborn
- A comfortable bedding
- An easy weighing
- A reliable measure
- Easy to clean



Weighing reliability and comfort

The Carewave by Médipréma cushion for baby-scales lays on the plate to welcome the new-born in good conditions during the weighing. The cold and hard contact of the plate which causes agitation of the baby is now replaced by a soft cocoon. The baby feels secure so it is quiet, and the measure of the weigh is easier and reliable. *Time saving for the care team and comfort for the new-born!*



Easy to position and secure

The Carewave by Médipréma cushion for baby-scales is light, easy to position. Made of polystyrene micro-granules, the ultra-flexible cushion stalls the baby without compressing him. Easy to clean, it will adapt to all types of baby-scales.

Specifications :

Dimensions	55 x 30 cm
Material inside the cushion	Micro-Granules in polystyrene
Matière enveloppe	Bi-elastic polyurethane slipcover
Cleaning	by sponging Water resistant Double seam
Packing	Per unit

A Pharmaouest product, distributed by: