



Night & Day canopy cover

for incubators Satis 3552 et 3600

Ref : 3658

For new-born comfort

- > Protect from surrounding light
- > Decrease noise pollution
- > Adapted to the environment of use

Decrease surrounding disturbance

Médipréma Night & Day canopy cover allows to reduce infant aggression feelings in a resuscitation unit. The new-born is protected from direct lighting and noise pollution is attenuated.



Adapted to the environment of use

Médipréma Night & Day canopy cover is very easy to place on an incubator and equally easy to pull back.

Inclined side can be turned up to keep partial visibility of the new-born. The cover has been designed in order not to stand in the way of sensors and pipes.

- > Complete or partial openings on all 4 sides
- > Opening on the head side to facilitate passage of tubes



Description :

Fabric	100% polyester (PES) Easy cleaning – Machine washable at 40°
Packaging	Per piece